BAUSCH+STROEBEL™

AFV 1010/2010

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AFV 1010/2010

FILLING AND CLOSING MACHINES FOR
AMPOULES AND VIALS

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The machines of this series are fully automatic filling and closing machines for both ampoules and vials. They are rotary machines processing 1 or 2 containers at each machine cycle and reaching an output of up to 3000 resp. 6000 containers per hour.

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* Supervisory Control And Data Acquisition

Multifunctional applications

Due to the compact design and accordingly small footprint, these machines can be found both in laboratories and in medium output production situations. With additional equipment they can be used alternatively for ampoules and for vials.

Short set-up time – higher production

Size changeover is simple and quickly carried out, so that these versatile machines can be used for different kinds of container and product.

Pharmaceutical design

Apart from GMP and FDA guidelines, standards applicable to the processing of pharmaceutical products provide the foundation for the design and operation of this machine. The arrangement of the work stations and the design above the work surface enhance a laminar air flow. Virtually all parts are made of corrosion-proof materials such as stainless steel and high-grade plastics. The design and surface finish ensure that the machine is easy to clean and care for.

In-line or stand-alone

Both machines can be equipped for use as stand-alone machines or for integration into a production line.

Top-range PLC

A programmable logic control (PLC) monitors and controls the machine functions on the basis of "Best quality at maximum output". The machine can be linked with data acquisition systems for the collection, evaluation and storage of data relevant to production (SCADA*).

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The advantages


- Universal application for both ampoules and vials
- Streamlined design for LF
- Gentle container transport
- Rotation during filling
- Precise dosing using valveless rotary piston pumps
- Sealing system with conventional gas burners or (for AFV 2010) with the innovative Laser Seal System
- Quick and easy size change

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Clearly arranged, versatile work area

- Loading containers from trays or automatic transfer from sterilizing tunnel upstream
- Intermittent transfer to the star wheel
- Container registration by photoelectric sensors to control following work steps (no container – no dosing – no stopper)
- Opening ampoules (if closed ampoules are processed)
- Nitrogen flushing before filling
- Precise dosing with valveless rotary piston pumps
- Nitrogen flushing after filling
- For ampoules: Nitrogen blanket during sealing
- Sealing by melting the stem and pulling off the ampoule tip
- For vials: Closing by inserting rubber stoppers
- Discharge into trays or transfer to a machine linked downstream

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Wider processing range with domed closed ampoules

Closed ampoules can also be processed if an opening station is added. This optional feature opens the ampoules directly prior to filling.

The first step is to melt a hole in the dome of the ampoule to equalize the pressure inside and out.

Then the hole is widened while the ampoule is rotated so that the aperture becomes circular.

Any condensate forming inside the ampoule stem is dried by the burners.

This two-step procedure ensures that the opening in the ampoule is ideal for filling the product and that the interior of the ampoule is in no way impaired.

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Valveless rotary piston pump(s) – self-priming and precise

The dosing station is equipped with either one or two valveless rotary piston pumps. These high-quality, precision pumps comprising only two parts guarantee utmost dosing accuracy. These pumps are simple to remove, install, clean and sterilize as are all the parts in contact with the product.

Depending on the product to be filled, either stainless steel or ceramic pumps can be supplied.

Dose adjustments within the pump volume can be made by increasing or reducing the pump stroke.

Both the central dose adjustment and the fine adjustment for each pump are made mechanically with setting spindles. They are equipped with digital position displays enabling the dosing quantities to be repeated, so that they can be reproduced with the same accuracy at any time.

An optional drip retraction device prevents problematic products from dripping at the end of the filling procedure.

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„Cleaning in place – sterilizing in place“ is a time-saving way of cleaning and sterilizing the product supply parts without having to remove them from the machine. All the parts in contact with the product are filled with the appropriate cleaning media and thoroughly rinsed.

Then they are sterilized with ultra-pure steam. When sterilization has been completed, the whole system is cooled and then the pressure in the system is maintained until the next production start.

The dosing needles are inserted in the CIP/SIP receptacle. All further position changes required for cleaning and sterilizing are either carried out manually after a prompt or are program controlled.

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
Tried and tested

The ampoules are sealed by fusing the glass and pulling off the tips of the ampoule stems.

The fuel gas and oxygen mixture produces a constant burner flame. The gas supply can be set individually on the respective flow meter for each work station.

If the flame goes out, the gas supply is shut off automatically.

The ampoule stems are heated in two steps and then as the glass is fused to a dome the tip is pulled off and discarded into a bin. The heat generated by the sealing process is drawn out of the work area by an extraction system.

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Sophisticated and innovative sealing system

With this sealing system a laser produces the energy to fuse the ampoule stems. The system produces no particles and requires few components above the work surface thereby guaranteeing optimal, unimpeded laminar air flow.

The ampoule stems are sealed in two steps. At the first station the fusing point is pre-heated and at the second station the tips are pulled off and discarded when the glass fusion temperature is reached. They drop into a collection bin.

The laser beam is only active while an ampoule is being sealed. During set-up work and if the transport position is vacant, the laser remains inactive. The laser module is largely maintenance-free.

The machine guard and the safety installations prevent the laser beams from escaping. Machine operation is in accordance with the laser hazard class 1 (according to DIN EN 60825-1).

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Upgrading to a universal machine for both ampoules and vials

With additional equipment for the supply and insertion of rubber stoppers both ampoules and vials can be filled and closed on the same machine.

The closures are fed to the closing station by a sorting vibrator and a track vibrator.

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Simple, quick and exact

Changing over the machine for different containers is quick and simple as complete size part sets are exchanged. The machine can then be adjusted to any container height.

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Ampoule processing using
laser technology

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Enhanced production reliability and greater versatility

- Special arrangement of the infeed and discharge areas
- Station for opening domed closed ampoules
- Automatic burner ignition
- Gas-flushing before and after filling and during sealing
- CIP/SIP equipment
- Centering device for ampoule stems
- Drip retraction
- Mobile version
- Program-controlled dose adjustment
- Various dosing vessels
- Discharge for IPC
- Sealing ampoules by laser (AFV 2010 only)
- Equipment to process vials
- LF unit
- Barrier technology
- Isolator
- A range of checking and alarm devices and counters
- Linkage to data acquisition systems for all data relevant to production (SCADA)
- Qualification (IQ/OQ)

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Subject to technical modifications.
Please note that the illustrations
may vary from the standard
version in some details.

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BAUSCH+STRÖBEL Maschinenfabrik Ilshofen GmbH + Co. KG
BAUSCH+STROEBEL Machine Company, Inc.