



Via XXV Aprile, 8 21054 Fagnano Olona (VA) - Italia Tel. +39 0331 1693557 email: inti@intisrl.it

Tube Filler Machines/Suppositories fillers Groninger MA200 De-nester

Images





Product details

Category:Tube Filler Machines/Suppositories fillersMachine:MA200 De-nesterMachine code:23-795

Manufacturer: Groninger

Year of construction: 1996

Description

Groninger MA200 De-nester for syringes from 1996, manual feeding of the nests filled with syringes into the nest conveyor (two rows side by side) and automatic feeding to the syringe outfeed station (two axes / horizontal X and vertical Y), where the syringes are removed from the nest in rows and pushed onto the syringe discharge conveyor (with the finger rests on the conveyor belt), the empty nests are conveyed out of the machine via a discharge chute and collected in a container, currently set up for syringes with a diameter of 10. 85 mm, output up to 400 syringes per minute, Siemens S7 PLC control with Siemens MP 377 12" touch panel, servo motors controlled, mains voltage 230/400 V(AC) / 50 Hz, electrical power 3 kW, compressed air supply network 7 bar